



OXIDATIVE STRESS, MEASUREMENT AND EVALUATION GUIDE

RELASSAY DIAGNOSTICS

2012



OXIDATIVE STRESS, MEASUREMENT AND EVALUATION

Free radicals are metabolic by products which are formed during energy conversion of nutrients in our bodies by using oxygen. Free radicals (reactive oxygen species) are in unstable form and in order to become stable they attack and damage cells. Antioxidants neutralize free radicals to protect cells from damage. Free radicals and antioxidants must be in balance thus antioxidants neutralize free radicals. If free radical level increased against antioxidant level oxidative damage occurs and this state is called oxidative stress.

Oxidative stress is not a disease, but it is a factor that can cause or accelerate the disease. Generally it is an important stimulus for preventive health measures. But the dangerous is oxidative stress is absence of any symptoms. If this condition is not diagnosed and treated, could result serious health problems and diseases. Exposure to toxins or pathogens, poor antioxidant defense system, irregular life style, extreme intense exercise and daily metabolic products, causes oxidative stress.

In some cases oxidative stress increase higher than normal levels.

For example;

Pregnancy

Hormone replacement therapy

Birth control pills

After heavy exercise

Excessive exposure to sunlight

Excessive alcohol consumption

Electromagnetic radiation

Smoking

Air Pollution

Malnutrition

Chronic inflammation

In such cases, must pay more attention to the antioxidant intake.

SERBEST RADİKALLERİN NEDEN OLDUĞU HASTALIKLARIN BAZILARI

GÖZ

- Diabetik retinopati
- Premature Retinopatisi
- Katarakt
- Yaşa Bağlı maküler Dejenerasyon

SOLUNUM YOLU

- Bronşiyal Astım
- Solunum Hasarı
- Yetişkin Solunum güçlüğü Sendromu

BATIN

- Akut gastrik Mukozal Lezyon
- Gastrik Ülser
- İskemik Enterit
- İskemik Kolit
- Karaciğer yağlanması

GENEL

- Diyabet
- Alerji
- Romatoid Hastalıklar
- Kanser
- Ateroskleroz

SEFALİK KISIM

- Parkinson
- Alzheimer
- Sinir Sistemi Hastalığı
- Travmatik Epilepsi

YÜZ


- Yaşlılık Lekesi
- Yaşlanma
- Kırışıklık
- Sarkma
- Güneş Lekesi

GÖĞÜS

- İskemik Aritmi
- Enfarktüs
- Hipertansiyon

ALT BATIN

- *Böbrek Yetmezliği





Reducing free radicals and increasing antioxidants by diet and supplements

There are many factors affecting the antioxidant defense system of the body. Individual genetic structure and environmental factors that body exposed. Unfortunately, the modern way of life, environmental pollution, stress, low-quality foods, unbalanced nutrition, causes remain to expose to more free radicals. So physiologically produced antioxidants are not enough to neutralize free radicals.

By rich in antioxidants diet defense of the body against free radicals can be increased.

Diet contains fruits, vegetables and oily seeds vitamins, antioxidants and other antioxidant micronutrients creates a good source of exogenous character and increases the cellular response to oxidative stress.


For example, essential vitamins such as vitamin C has suppressive and reducer effect on many different reactive oxygen species and not synthesized by human body. Epidemiological studies support that people who eat rich in antioxidants healths are better. The World Health Organization (WHO) recommends as a result of epidemiological studies, daily consumption of at least 400g of fruit and vegetables.

Every part of the cell is maintained by a different antioxidants. Therefore fed on by a variety of nutrition is important. For example some of the food items protects against free radicals inside the cell, outside the cell wall, inside the blood surrounding the cell and the cell membrane. For this reason, we need both exogenous and endogenous antioxidants

As examples of antioxidants:

- Vitamin C, protects vitamin E and refreshes the extracellular fluid.
- Vitamin E protects cell wall.
- Carotenoids (In the content of the fruits and vegetables that people consumed, founded over 60 caretenoids, beta-carotenoids, laykopen, lutein and other colored compounds, yellow, red, orange.)
- Thiols protects inside of the cell (sulfur-containing compounds such as lipoic acid, and glutathione) .
- Co-enzyme Q10 protects mitochondria and slows aging.
- Flavonoids protects DNA, elastin and collagen tissue, thus slow down the aging process.
- Enzymes including superoxide dismutase, glutathione peroxidase, catalase, can neutralize free radicals in the cells. They can not gotten from the outside they are endogenous antioxidants produced by the body. However, they need the dietary selenium, copper, zinc and trace minerals such as manganese.

Making frequent changes in the diet may not enough to organize antioxidant intake. Diet should be prepared by an expert and necessary foods contains antioxidants should added to the diet.



General recommendations to reduce oxidative stress

- Get a balanced diet, consume adequate amounts of fruits and vegetables, if possible, be sure to choose organic of them. Drink plenty of water, avoid fatty, fried and processed foods.
- Avoid excessive consumption of sugar and salt, and choose products that are not refined.
- Limit your alcohol intake.
- Consume more tea (green tea), and less coffee.
- Do mid-level exercise on a regular basis, but do not overdo it. Much forcing your body, increase the level of free radicals.
- Stay as far away as possible from pollution, toxic materials causes oxidative stress. Stay away from industrial pollution, car pollution, cigarette smoke, and toxic substances in food.
- Insomnia causes oxidative stress, thus try to get regular sleep.
- Avoid long-term exposure to ultra-violet rays and being under the sun. Apply strong protection creams under the sun.
- Regularly measure your blood sugar, cholesterol and blood pressure and keep under control.
- If you are taking hormones (birth control pills or hormone replacement therapy) keep control of your oxidative stress.
- Keep away from excessive mental and physical stress.
- Give yourself time to relax and do things that make you happy.

**OXIDATIVE STRESS MEASUREMENT**

Serum (or plasma) concentrations of different antioxidants can be measured in laboratories separately, but the measurements are time-consuming, labor-intensive, costly, and require complicated techniques. Because the measurement of different antioxidant molecules separately is not practical and their antioxidant effects are not additive.

With Rel Assay Total Antioxidant Status(TAS) and Total Oxidant Status (TOS) kits you can easily measure non-enzymatic all antioxidant and oxidant molecules.

For measurement procedures and booklets you can visit www.rel assay.com, send e-mail to rel assay@rel assay.com or make a phone call to +90 342 220 49 87.

Collection of samples

For this measurement serum is used as sample. Hemolyzed or lipemic serum is not required. The sample should be taken after fasting for 10-12 hours. Sample must not taken after exercise.

Storage of samples

Serum 2-8 ° C 24 hours

Serum -20 ° C 3 months

Serum -80 ° C 1 Year

TAS Test Principle

Antioxidants in the sample reduce dark blue-green colored ABTS radical to colorless reduced ABTS form. The change of absorbance at 660 nm is related with total antioxidant level of the sample. The assay is calibrated with a stable antioxidant standard solution which is traditionally named as Trolox Equivalent that is a vitamin E analog .

TOS Test Principle

Oxidants present in the sample oxidize the ferrous ion–chelator complex to ferric ion. The oxidation reaction is prolonged by enhancer molecules, which are abundantly present in the reaction medium. The ferric ion makes a colored complex with chromogen in an acidic medium. The color intensity, which can be measured spectrophotometrically, is related to the total amount of oxidant molecules present in the sample. The assay is calibrated with hydrogen peroxide and the results are expressed in terms of micromolar hydrogen peroxide equivalent per liter (μmol H₂O₂Equiv./L)





Evaluation of the results

Rel Assay TAS and TOS test results can be evaluated independently in accordance with the following tables or we recommend to put the results of both to 9-diagnosquare to evaluate.

TAS REFERANS DEĞERLER		
(mmol Trolox Equiv./L)		
>2,0		Çok iyi
1,45	2,00	Normal
1,20	1,45	Normal Kabul Edilebilir
1,00	1,20	Düşük Antioksidan Seviyesi
<1,20		Çok Düşük Antioksidan Seviyesi

TOS REFERANS DEĞERLER		
(μmol H ₂ O ₂ Equiv./L)		
<5,00		Çok iyi
8,00	5,00	Normal Değer
12,00	8,00	Yüksek Oksidan Seviyesi
>12,00		Çok Yüksek Oksidan Seviyesi

- This reference range is prepared in the light of the scientific data in accordance with the international literature.
- Samples used in this study carried out in Turkey by the Turkish population.
- In this study, samples were taken from 20-60 years of age healthy individuals has taken place as "normal" status.
- It is recommended to evaluate the results of TAS and TOS according to the 9-DiagnoSquare .
- It should be noted that each laboratory determine its own reference range. This study is a sample, prepared by the manufacturer by itself. The manufacturer does not accept responsibility for any problems arising from the use of this data.



1 - VERY HIGH Oxidative Stress:

2 - HIGH OXIDATIVE STRESS:

3 - OXIDATIVE STRESS:

4 - IMMUNITY AND METABOLISM DISORDERS:

5 - BALANCE – MAY BE THE BEGINNING OF A DISEASE

6 - NORMAL LEVEL


7 - GOOD LEVEL

8 - VERY GOOD LEVEL

9 - BEST LEVEL




Oxidative Stress




1 - VERY HIGH Oxidative Stress:

Oxidants too much, antioxidants is very low, this is quite dangerous. Individuals in this region must be suspected from degenerative diseases, cardiovascular diseases and cancer diseases and advanced examinations must be consulted. In addition, use of medication, chemotherapy and radiation therapy can cause very high oxidative damage. The individual must be examined in detail, besides specifically treatment must be supported by antioxidants. TAS and TOS levels must be checked every week.




2 - HIGH OXIDATIVE STRESS:

The state of high oxidant level beside normal antioxidant level. This may be initial phase of many degenerative diseases. Hormone replacement therapy, birth control pills, sedentary lifestyle, alcohol, smoking and even obesity can cause high-oxidative level. In addition it can be seen in pregnancy. In addition, high sporting activity can cause this case. Patient history must be listened carefully and, if necessary supported by other tests. Oxidant reducer treatment must be executed and TAS and TOS levels must be checked for every 15 days.




3 - OXIDATIVE STRESS:

The condition with normal levels of oxidants, but with not enough antioxidants. In these individuals due to malnutrition, attention deficit, infectious diseases, the slow pace of growth and development in children, mental retardation, premature aging and such as diseases can be seen. Individuals' lifestyles and eating habits could be examined in detail and they should be supported with antioxidants by a dietitian who has results in this area. TAS and TOS tests should be repeated after 15 days.




4 - IMMUNITY AND METABOLISM DISORDERS:

Oxidants and antioxidants are both low than normal, this situation is rare. Decrease in the level of oxidants can be interpreted good but decrease in antioxidants is not a good situation at all. The results can be like this at the breakdowns of immune system. Metabolic diseases should be considered in such situations, and further investigations should be requested. TAS and TOS tests should be repeated after 15 days.




5 - BALANCE – MAY BE THE BEGINNING OF A DISEASE


Oxidants and antioxidants both high, this situation is rare. This occurs usually after intense exercise in athletes. This condition may be normal for athletes but for healthy individuals increase of oxidants that much can cause some problems. In healthy individuals this can be observed at pre-disease cases or cases of acute exposure to poor environmental conditions. Required specific therapy must be implemented and to confirm decrease in oxidant level, test should be repeated after 15 days.

**6 - NORMAL LEVEL**


Previous studies have shown that, under normal conditions Oxidants and Antioxidants is in equilibrium for the organism, and the vast majority of healthy individuals are in this group. Nowadays, enough information is available about the importance of antioxidant defense system and the hazards of oxidants. Therefore should pay attention to lifestyle to not to slide this balance toward oxidative stress and test should be repeated every 2 months.

**7 - GOOD LEVEL**

The normal level of oxidants, and higher antioxidant level, in this situation there is no oxidative stress. Generally individuals who care about diversity of plant foods in their diet have this results. In addition, this situation can be found for the individuals who takes antioxidant support.

**8 - VERY GOOD LEVEL**

Oxidants is at the below of normal level, antioxidants is in normal level, in this situation there is no oxidative stress. Decrease in oxidants could be related to lifestyle and living environment. Generally, found in individuals who is far away from city life and naturally nourished.

**9 - BEST LEVEL**

Above normal level of antioxidants, below normal level of oxidants, this situation is described as the best level. If the individual is not taking antioxidant supplements and reached to this situation means individual is at disease risk-free group.

MEGA TIP SANAYİ VE TİCARET LİMİTED ŞİRKETİ

BAHÇELİEVLER MAHALLESİ KAYMAKAM İSMAILPAŞA SOKAK NO:10 ŞAHİNBEY / GAZİANTEP / TÜRKİYE

TELEFON: +90 342 230 37 52 FAKS: +90 342 220 49 87

WWW.RELASSAY.COM – RELASSAY@RELASSAY.COM



References;

- Halliwell B, Gutteridge JMC. Free radicals in biology and medicine. Third ed. Oxford: Oxford Science Publications; 2000. p. 617–24.
- Tarpey MM, Wink DA, Grisham MB. Methods for detection of reactive metabolites of oxygen and nitrogen: in vitro and in vivo considerations. *Am J Physiol, Regul Integr Comp Physiol* 2004;286(3):R431–44.
- Yeni E, Gulum M, Selek S, et al. Comparison of oxidative/antioxidative status of penil corpus cavernosum blood and peripheral venous blood. *Int J Impot Res* 2005;17(1):19–22.
- Yanik M, Erel O, Kati M. The relationship between potency of oxidative stress and severity of depression. *Acta Neuropsychiatr* 2004;16(4):200–3.
- Harma M, Harma M, Erel O. Increased oxidative stress in patients with hydatiform mole. *Swiss Med Wkly* 2003;133(41–42):563–6.
- Nakamura K, Endo H, Kashiwazaki S. Serum oxidation activities and rheumatoid arthritis. *Int J Tissue React* 1997;9(4):307–16.
- Ceylan E, Gulsun A, Gencer M, Aksoy N. A new parameter in the detection of tuberculosis activity: reactive oxygen metabolites. *Respiration* 2005;72(2):156–9.
- Erel O. A new automated colorimetric method for measuring total oxidant status. *Clin Biochem.* 2005 Dec;38(12):1103–11.
- Dyt. Andaç Yeşilyurt
- Rel Assay Research Laboratory