PARAOXONASE

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10.KAYNAKÇA

1. Ekmekci A. Koroner arter hastalığında ve tedavisinde damar endotel fonksiyonları. **CIBA 1995**, İstanbul.
2. Onat A. Risk factors and cardiovascular disease in Turkey. *Atherosclerosis* 2001; 156:1-10
3. Kültürsoy H. Koroner kalp hastalığı primer ve sekonder korunma. **ARGOS 2001**; T.A.Ş İstanbul
4. Koylan N. Koroner arter hastalığı epidemiyolojisi lipid düşürücü tedavi ile ilgili büyük klinik çalışmalar. *Türk Kardiyoil Dern Arşv* 2000;13:9-20
5. Castelli W P, Wilson P W F, Deny D, Anderson K D. Cardiovascular risk factors in the elderly. *Am J Cardiol.* 1989; 63:12-19
6. Gordan T, Castelli W P, Hjortland M C. Lipoproteins cardiovascular diseases and death. The Framingam Study. *Arch Intem Mad.* 1981; 141: 1128
7. Mahley R W, Innearity T L, Rall S.C, Weisgraber K.H, Plasma lipoproteins, apolipoprotein structure and function. *Journal of Lipid Research.* 1984; 25: 1277-94
8. Mackness M I, Arrol S, Durrington P N, Paraonase prevents accumulation of lipoperoxides in low-density lipoprotein. *FEBS Lett.* 1991; 286: 152-154,
9. Hernandez A F, Pla A, Valenzuela A., Gil F, Hougen H P, Villanueva E. Paraonase activity in human pericardial fluid its relationship to coronary artery disease. *Int J Leg Med.* 1993; 105: 321-324
10. Kelso G J, Stuart W D, Richter R J, Furlong C E, Jordan-Starck T C, Harmony J A. Apolipoprotein J is associated with paraonase in human plasma. *Biochemistry* 1994; 33:832-839
11. Mackness M I, Arrol S, Abbott C, Durrington P N. Protection of low density lipoprotein against oxidative modification by high-density lipoprotein associated paraonase. *Atherosclerosis* 1993; 104:129-135
12. Mackness M I, Harty D, Bhatnagar D, Winocour P R, Flirli S, Ishola M, Durrington P. Serum paraonase activity in familial hypercholesterolemia and insulin-dependent diabetes mellitus. *Atherosclerosis* 1991; 86:193-199
13. Eckerson H W, Wyte C M, La Du B N. The human serum paraonase/arylesterase polymorphism. *Am J Hum Genet.* 1983; 35:1126-1138
14. Navab M, Hama-Levy S, Van Lenten R I, Fonarow G C, Cardinez C J, Castellani L W, Brennan M L, Lusis A I, Fogelman A M, La Du B N. Mildly oxidized LDL induces an increased apolipoprotein J/paraonase ratio. *J Clin Invest.* 1997; 99:2005-2019
15. Odowara M, Tachi Y, Yamashita K. PON1 polymorphism (Gln 192-Arg) is associated with coronary heart disease in Japanese non-insulin-dependent diabetes mellitus. *J Clin end Metb* 1997; 82:2257-2260
16. Heinecke I W, Lusis A J. Paraonase-gene polymorphisms associated with coronary heart disease: support for the oxidative damage hypothesis? *Am J Hum Genet* 1998; 62:20-24
17. Kao Y L, Donaghue K, Chan A, Knight J, Silink M. A variant of paraonase (PON1) gene is associated with diabetic retinopathy in IDDM. *J Clin Endocrinol Metab* 1998; 83:2589-2592
18. Gan K N, Smolen A, Eckerson H W, La Du B N. Purification of human serum paraonase/arylesterase. Evidence for one esterase catalyzing both activities. *Drug Metab Dispos.* 1991; 19:100-106
19. Costa L G, Mc Donald B E, Murphy S D, Omenn G S, Richter R J. Serum PON and its influence on Paraon. *Tox and Applied Pharmaeology* 1990; 103:66-76
20. Hegele R A, Connelly P W, Scherer S W, Hanley A J, Harris S R, Tsui L C, Zinman B. Paraonase 2 gene (PON2) G148 variant associated with elevated fasting plasma glucose in noninsulin-dependent diabetes mellitus. *J Clin Endocrinol Metab* 1997; 82:3373-3377
21. Costa L, Vitalone A, Cole T B, Furlong C E. Modulation paraonase (PON1) activity. *Biochemical Pharmacology* 2004
22. James R, Leviev I, Ruiz J, Passa P, Fuegel P, Gavin MCB. Promoter polymorphism T(-107)C of the paraonase PON1 gene is a risk factor for coronary heart disease in type 2 diabetic patients. *Diabetes* 2000; 29:390-39
23. Paragh G, Balla P, Katona E, Seres I, Egerhazi A, Degrell I, Serum paraonase activity changes in patients with Alzheimer's disease and vascular dementia. *Eur Arch Psychiatry Clin Neurosci* 2002; 252:63-67
24. Ferre N, Marsilach J, Camps J, Rull A, Coll B, Tous M, Joven J. Genetic association of paraonase-1 polymorphisms and chronic hepatitis C virus infection. *Clinica chimica acta.* 2005; 18:112-118
25. Min J, Park H, Park B, Kim YJ, Park J, Lee H, Ha E, Park E A, Hong Y C. Paraonase gene polymorphism and vitamin levels during pregnancy: relationship with maternal oxidative stres and neonatal birtweights. *Reproductive toxicology* 2006; 147:363-367
26. Clarimon J, Eerola J, Hellstrom O, Tienari PJ, Singleton A, Paraonase 1 (PON1) gene polymorphisms and Parkinson's disease in a Finnish population. *Neurosci Lett* 2004; 367:168-170
27. Paşca S.P, Nemeş B, Vlase L, Gaygi C.E, Dronca E, Miu A.C, Dronca M., High levels of homocysteine and low serum paraonase 1 aylesterase activity in chidren with autism. *Life sciences* 2006; 78:2244-2248
28. Raiszadeh F, Solati M, Etemadi A, Azizi F. Serum paraonase activity before and after treatment of thyrotoxicosis. *Clin Endocr* 2004; 60:75-80
29. Delgado Alves J, Ames PR, Donohue S, Stanyer L, Nourooz-Zadeh J, Ravirajan C, Isenberg DA, Noorouz-Zadeh J. Antibodies to highdensity lipoprotein and β -2 glycoprotein-I are inversely correlated with paraonase activity in systemic lupus erythematosus and primary antiphospholipid syndrome. *Arthritis Rheum* 2003; 48:284
30. Mackness MI, Mackness B, Arrol S, Wood G, Bhatnagar D, Durrington PN. Presence of paraonase in human interstitial fluid. *FEBS Letters* 1997; 416:377-380.
31. Griffith MK, Virella GT, Stevenson HC, Lopes-Virella MF. Low density lipoprotein metabolism by human macrophages activated with low density lipoprotein immune complexes. A possible mechanism of foam cell formation. *J Exp Med* 1988; 168:1041-1059
32. Szafran Z, Nowak J, Szafran H, Janik A. Esterolytic activity of blood serum in infants with hypertrophic pylonic stenosis. *J Clin Chem Clin Biochem* 1997; 17:321-324
33. Van Lenten BJ, Hama SY, de Beer FC, Staffonini DM, McIntyre TM, Prescott SM, La Du BN, Fogelman AM, Navab M. Anti-inflammatory HDL becomes proinflammatory during the acute phase response. *J Clin Invest* 1995; 96:2758-2767

34. Berliner JA, Navab M, Fogelman AM. Atherosclerosis: basic mechanisms. Oxidation, inflammation, and genetics. *Circulation* 1995; 91:2488–2496
35. Marchesani M, Hakkarainen A, Tuomainen T P, Kaikkonen J, Pukkala E, Uimari P, Seppalo E, Matikainen M, Kallioniemi O P, Schleutker J, Lehtimäki T, Salonen J T. New paraoxonase 1 polymorphism I102V and the risk of prostate cancer in Finnish men. *J Natl Cancer Inst* 2003; 95:812–818
36. Erel O. A new automated colorimetric method for measuring total oxidant status. *Clin Biochem* 2005; 38:1103–1111
37. Rel Assay Paraoxonase Assay Kit (MEGA TIP SAN VE TIC LTD STI) <http://www.relassay.com/?sf=3&M=4&IDetay35=38>
38. SELEK S. Harran University Faculty of Medicine, Sanliurfa, Turkey.
39. Mackness M, Durrington P, Mackness B. Paraoxonase 1 activity, concentration and genotype in cardiovascular disease. *Curr Opin Lipidol* 2004; 15:399–404.
40. suchocka Z, swatowska J, pachecka J, suchocka P; RP-HPLC determination of Paraoxonase activity in human blood serum; *J of pharmaceutical and biomedical analysis* 2006; 42:113–119
41. Lipincott W. Paraoxonase a cardioprotective enzyme: continuing issues, *Curr Opin Lipidol* 2004; 15:261–267
42. Shamir R, Hartman C, Karry R, Pavlotzky E, Eliakim R, Lachter J, Suissa A, Aviram M; Paraoxonases 1, 2 and 3 are expressed in human and mouse gastrointestinal tract and in Caco-2 cell line: selective secretion of PON1 and PON2; *Free Radical biology* 2005; 39:336–344
43. Yang X, Gao Y, Zhou J, Zhen Y, Yang Y, Wang J, Song L, Liu Y, Xu H, Chen Z, Hui R; Plasma homocysteine thiolactone adducts associated with risk of coronary heart disease; *C Chimica acta* 2006; 364: 230–234
44. Gülcü F, Gürsu F; The standardization of paraoxonase and arylesterase Activity Measurements; *Turkish J Biochem* 2003; 28:45–49
45. Aviram M, Davies J. A, Paraoxonase 1 and atherosclerosis: Is the gene or the protein more important? *Free Radical Biology* 2004; 37:1317–1323
46. La Du BN, Eckerson HW. The polymorphic paraoxonase/arylesterase isozymes of human serum. *Fed Proc* 1984; 43:2338–2341
47. Karakaya A., Suzen S., Sardas, S, Karakaya AE, Vural N. Analysis of the serum paraoxonase/arylesterase polymorphism in a Turkish population. *Pharmacogenetics* 1991; 1:58–61
48. Mallinckrodt M, Geldmacher V, Hommel G, Dumbach J. On the genetics of the human serum paraoxonase. *Hum Genet* 1979; 50:313–326
49. Kleemola P, Freese R, Jauhiainen M, Pahlman R, Altfhan G, Mutanen M; Dietary determinants of serum paraoxonase activity in healthy humans. *Atherosclerosis* 2002; 160:425–432
50. Ombres D, Pannitteri G, Montali A, Candeloro A, Seccareccia F, Campagna F, Cantini R, Campa P P, Ricci G, Arca M. The Gln-Arg192 Polymorphism of Human Paraoxonase Gene Is Not Associated With Coronary Artery Disease in Italian Patients. *Atheroscler Thromb Vasc Biol.* 1998; 18:1611–1616
51. Garin MCB, James RW, Dussoix P, Blanchè H, Passa P, Froguel P, Ruiz J. Paraoxonase polymorphism Met-Leu 54 is associated with modified serum concentrations of the enzyme: a possible link between the paraoxonase gene and increased risk of cardiovascular disease in diabetes. *J Clin Invest* 1997; 99:62–64.
52. Agachan B, Ergen H A, Karaali Z E, Isbir T, PON1 55 and 192 Polymorphism and Its Effects to Oxidant-Antioxidant System in Turkish Patients with Type 2 Diabetes Mellitus. *Physiol Res* 2005; 54:287–293
53. Karakaya A, Ibis S, Kural T, Kose SK, Karakaya A.E, Serum paraoxonase activity and phenotype distribution in Turkish subjects with coronary heart disease and its relationship to serum lipids and lipoproteins. *Chem Biol Interact* 1999; 118:193–200
54. Furlong C E, Li W F, Brophy VH, Jarvik G P, Richter R J, Shih D M, Lusis A L, Costa L G: The PON1 gene and detoxication. *Neurotoxicology* 2000; 21:581–587
55. Furlong CE, Richter RJ., Seidel S L, Motulsky AG. Role of genetic polymorphism of human plasma paraoxonase/arylesterase in hydrolysis of the insecticide metabolites chlorpyrifos oxon and paraoxon. *Am J Hum Genet* 1998; 43:230–238
56. Mackness M I, Arrol S, Durrington P N, Paraoxonase prevents accumulation of lipoperoxides in low-density lipoprotein. *FEBS Lett.* 1991; 286: 152–154,