

TAS REFERENCE VALUES

(mmol Trolox Equiv./L)

>2,0		Very good
1,45	2,00	Normal
1,20	1,45	Tolerable
1,00	1,20	Low Antioxidant Level
<1,20		Very Low Antioxidant Level

- This reference range study was prepared in consideration of scientific data in accordance with the international literature.
- This study carried out in Turkey with samples of Turkish population.
- In this study, samples were taken from healthy individuals of 20-60 years of age has taken place at "normal" status.
- TAS results alongside TOS results should be analyzed and evaluated according to the 9-DiagnoSquare.
- It should be noted that each laboratory should determine its own reference range.
- This study is just a sample that prepared by the manufacturer for itself.
- The manufacturer does not accept any responsibility for problems resulting from the use of this data.